

medipal

# maceratable dry wipes



High absorbency  
Less product usage  
Less waste



Made from  
100% natural materials  
Eco-friendly  
FSC - certified fibres



Maceratable  
Biodegradable  
100% recyclable packaging

## Products in this range



**Dimensions** – 280mm x 250mm per wipe

**Count** – 75 wipes per pack

**Primary packaging** – Flow wrap 100% recyclable CPP film

**Packs per case** – 16

**Cases per pallet** – 36

**Product Code** – DW600110MP

pal®

sales@palinternational.com | +44(0)1455 555 700 | palinternational.com



## Gentle, versatile & environmentally friendly

Medipal Maceratable Dry Wipes are designed for everyday cleansing and refreshing. Soft and gentle, they are ideal for fragile and sensitive skin types, offering total peace of mind in healthcare settings.

- ✔ **Soft & gentle** – Specially formulated for sensitive skin
- ✔ **Highly absorbent & large size** – Ideal for full-body cleansing
- ✔ **Maceratable** – Suitable for maceration

## Disposal

Discard used wipes in the appropriate waste bin following your local guidelines. Suitable for maceration. Biodegradable.

## How to use

### Step 1 Before you start

Pay particular attention to the elderly, children and patients with sensitive skin. Select your Medipal Patient Care Product. Use a fresh wipe for each step. Follow your standard moving and handling techniques when bed bathing patients. Do not use on wounds or broken skin. Please follow your local guidance for correct PPE usage.



### Step 2 8 Step bathing technique



**1. Head & neck**  
Use the wipe on face, paying particular attention to nose, ears and neck but avoid eyes.



**2. Left arm & armpit**  
Starting from upper shoulder, clean down first arm and hand. Pay particular attention to areas between the fingers, elbow joint, and any skin folds. Clean under armpit.



**3. Right arm & armpit**  
Starting from upper shoulder, clean down first arm and hand. Pay particular attention to areas between the fingers, elbow joint, and any skin folds. Clean under armpit.



**4. Left leg & foot**  
Starting from hip joint, clean down first leg to the foot and sole. Pay particular attention to areas between the toes, knee joint and any skin folds.



**5. Right leg & foot**  
Starting from hip joint, clean down first leg to the foot and sole. Pay particular attention to areas between the toes, knee joint and any skin folds.



**6. Chest & abdomen**  
Clean the chest being careful with the sensitive areas.



**7. Back & lower back**  
Clean the back and lower back.



**8. Groin & pubic area**  
Wash the groin and perineal area following your standard nursing practice.

### Step 3 After care Allow the patients skin to air dry.



**Pal International Limited**  
Unit 3 Mountpark,  
Bardon II, Victoria Road,  
Ellistown,  
Coalville LE67 1FA, UK  
E: info@palinternational.com

**EU Representative:**  
Pal Hygiene Products Limited  
Unit 5B & Unit 5H, Fingal Bay Business Park,  
Balbriggan, Co. Dublin.  
Republic of Ireland, K32 NY57  
E: info@palinternational.com