# medipal



# skin cleansing wash mitts & cloths







Lanolin free Paraben free Deodorising



Biodegradable wash cloths
Packaging 100% recyclable
Made in Britain



Skin friendly pH5.5 Refreshing

Hydrating



#### **Products in this range**





# **8 Skin Cleansing Wash Mitts,** Flow Wrap

Packs per case: 25

Product code: S658110MP

8 Skin Cleansing Wash Cloth,

Flow Wrap

Packs per case: 30

Product code: S658110WPFSB

#### Chemical formulation:

Medipal Skin Cleansing solution is designed by our in house experts to help keep skin clean, fresh and nourished. The formulation is maintained at a skin balanced pH and is dermatologically tested.

#### **Materials:**

Spunlace material is made by mechanically entangling fibres through the application of high pressure water jets. This gives the material a thick, soft and strong feel. The apertured texture gives the wipes a high liquid capacity which makes them suitable for wiping larger surfaces and the drape of the material gives a 'cloth-like' feel.







## Cleanse, hydrate, moisturise

Medipal Skin Cleansing Wash Cloths and Mitts provide gentle, effective cleansing while hydrating and soothing the skin, ideal for patient hygiene in healthcare settings.



- Alcohol-free & pH balanced
- Hydrating & moisturising

## How to use

Step 1

Before you start

Pay particular attention to the elderly, children and patients with sensitive skin. Select your Medipal Patient Care Product. Use a fresh wipe for each step. Follow your standard moving and handling techniques when bed bathing patients. Do not use on wounds or broken skin. Please follow your local guidance for correct PPE usage.



Step 2

### 8 Step bathing technique



1. Head & neck

Use the wipe on face, paying particular attention to nose, ears and neck but avoid eyes.



2. Left arm & armpit

Starting from upper shoulder, clean down first arm and hand. Pay particular attention to areas between the fingers, elbow joint, and any skin folds. Clean under armpit.



Right arm & armpit

Starting from upper shoulder, clean down first arm and hand. Pay particular attention to areas between the fingers, elbow joint, and any skin folds. Clean under armpit.



4. Left leg & foot

Starting from hip joint, clean down first leg to the foot and sole. Pay particular attention to areas between the toes, knee joint and any skin folds.



5. Right leg & foot

Starting from hip joint, clean down first leg to the foot and sole. Pay particular attention to areas between the toes, knee joint and any skin folds.



6. Chest & abdomen

Clean the chest being careful with the sensitive areas.



7. Back & lower back Clean the back and lower back.



8. Groin & pubic area

Wash the groin and perennial area following your standard nursing practice.



After care



**Disposal** 



Allow the patients skin to air dry.



macerate

flush
Do not



Discard used wipes in the appropriate waste bin following your local guidelines.



Unit 3 Mountpark, Bardon II, Victoria Road, Ellistown, Coalville LE67 1FA, UK E: info@palinternational.com

#### **EU Representative:**

Pal Hygiene Products Limited Unit 5B & Unit 5H, Fingal Bay Business Park, Balbriggan, Co. Dublin. Republic of Ireland, K32 NY57 E: info@palinternational.com



