



medipal

Disinfectant Wipes



12 products in this range

- Intended for disinfecting surfaces of non-invasive medical devices
- Dual-active chemistry for broad spectrum efficacy
- Dermatologically tested
- Intermediate level disinfectant



Disinfectant Wipes

Key Highlights

- Best selling range of everyday disinfection wipes
- Ready to use format provides fast, reliable disinfection in healthcare environments
- Alcohol & PHMB Free
- Dermatologically tested
- Available in flow wraps, cannisters and buckets
- 89% plastic reduction for refill bag vs the bucket with lid
- Tritex material, canister, buckets, lids, and labels are fully recyclable
- Solution is biodegradable
- Cardboard used in our outer packaging is FSC approved

How to use


STEP 1 Risk assessment

Please follow your agreed risk assessment policy guidelines regarding the use of PPE. Check medical device equipment manufacturers guidelines for DDAC and Triamine suitability.




STEP 2 Wipe selection

Check pack for use period and dispense a wipe. Close pack to protect remaining wipes. If visible soil is present, remove with a wipe. Use a fresh wipe to disinfect the surface.



STEP 3 S-Shape technique

Wipe the surface in an S-Shape moving from clean to dirty. Use the wipe flat not scrunched, ensuring the entire surface is wetted. Do not go over the same area twice with the same wipe. Use a fresh wipe if your wipe becomes soiled or dry.




STEP 4 Discard

Do not re-use wipes. Discard used wipes in the appropriate waste bin or macerator, following your local agreed guidelines.



STEP 5 Let dry

Allow the surface to dry naturally before use. Leave the surface for the stated contact time.



Use biocides safely. Always read the label and product information before use.

Technical details

Chemical Formulation

High specification formulation

Dual-active ingredients provide broad spectrum efficacy

Formulated with surface active agents

Substrate

Pal's proprietary Tritex® material is made by bonding layers of spunbond and meltblown polypropylene fibres together in a high heat process

The addition of a hydrophilic coating means optimum and consistent liquid discharge

Tritex® is 100% synthetic

Substrate key features



High Strength



Excellent solution retention



Very low linting



Efficient solution release

Not intended for end point sterilisation of medical devices

Efficacy details

Standard	Organism Group	Organism(s)	Contact Time	Conditions
EN 1276	Bactericidal	<i>Pseudomonas aeruginosa</i> <i>Staphylococcus aureus</i> <i>Enterococcus hirae</i> <i>Escherichia coli</i>	1 minute	Clean
EN 13727	Bactericidal	<i>Pseudomonas aeruginosa</i> <i>Staphylococcus aureus</i> <i>Escherichia coli</i> K12 <i>Enterococcus hirae</i>	1 minute	Clean
EN 13727	Bactericidal	<i>Pseudomonas aeruginosa</i> <i>Staphylococcus aureus</i> <i>Enterococcus hirae</i>	1 minute	Clean
EN 13727	Bactericidal	<i>Pseudomonas aeruginosa</i> <i>Staphylococcus aureus</i> <i>Enterococcus hirae</i>	1 minute	Dirty
EN 16615	Bactericidal	<i>Pseudomonas aeruginosa</i> <i>Staphylococcus aureus</i> <i>Enterococcus hirae</i>	1 minute	Clean
EN 16615	Bactericidal	<i>Pseudomonas aeruginosa</i> <i>Staphylococcus aureus</i> <i>Enterococcus hirae</i>	1 minute	Clean
EN 16615	Bactericidal	<i>Pseudomonas aeruginosa</i> <i>Staphylococcus aureus</i> <i>Enterococcus hirae</i>	1 minute	Dirty
EN 14348	Mycobactericidal	<i>Mycobacterium avium</i>	1 minute	Clean
EN 14348	Mycobactericidal	<i>Mycobacterium terrae</i> <i>Mycobacterium avium</i>	1 minute	Clean
EN 1650	Fungicidal Yeast (C.a) Fungi (A.b)	<i>Candida albicans</i> <i>Aspergillus brasiliensis</i>	1 minute	Clean
EN 13624	Yeasticidal	<i>Candida albicans</i>	1 minute	Clean
EN 13624	Fungicidal Yeast (C.a) Fungi (A.b)	<i>Candida albicans</i> <i>Aspergillus brasiliensis</i>	2 minutes	Clean
EN 13624	Yeasticidal	<i>Candida albicans</i>	1 minute	Dirty
EN 13624	Yeasticidal	<i>Candida auris</i>	1 minute	Clean

Efficacy details

Standard	Organism Group	Organism(s)	Contact Time	Conditions
EN 13697	Yeasticidal	Candida albicans	1 minute	Dirty
EN 16615	Yeasticidal	Candida albicans	1 minute	Clean
EN 16615	Yeasticidal	Candida albicans	1 minute	Clean
EN 16615	Yeasticidal	Candida albicans	1 minute	Dirty
EN 14476	Virucidal (Full)	Murine norovirus, strain s99 Berlin	1 minute	Clean
		Poliovirus type 1, strain LSc2ab	2 minutes	Clean
EN 14476	Virucidal (Human coronavirus)	Human coronavirus 229E	1 minute	Clean
EN 14476	Virucidal (Human coronavirus)	Human coronavirus 229E	15 & 30 seconds	Dirty
EN 14476	Virucidal (Limited)	Adenovirus type 5, strain Adenoid 75 Murine norovirus, strain s99 Berlin	1 minute	Clean
EN 14476	Virucidal (Enveloped)	Vaccinia virus VR-1508 (Modified Vaccinia Ankara)	1 minute	Dirty
EN 14476	Virucidal (Limited)	Murine norovirus, strain s99 Berlin	2 minutes	Dirty
		Adenovirus type 5, strain Adenoid 75	2 minutes	Dirty
EN 14476	Virucidal (Limited)	Adenovirus type 5, strain Adenoid 75 Murine norovirus, strain s99 Berlin	1 minute	Clean
EN 16777	Virucidal (Human coronavirus)	Human coronavirus 229E	1 minute	Clean
EN 16777	Virucidal (Limited)	Adenovirus type 5, strain Adenoid 75 Murine norovirus, strain s99 Berlin	1 minute	Clean

Additional Organisms

Standard	Organism Group	Organism(s)	Contact Time	Conditions
EN 13727	Bactericidal	<i>Acinetobacter baumannii</i> (Multidrug-resistant) <i>Enterococcus faecalis</i> (Vancomycin-resistant) <i>Enterococcus faecalis</i> (Vancomycin-resistant) <i>Enterococcus faecium</i> (Vancomycin-resistant) <i>Escherichia coli</i> (Carbapenem-resistant) <i>Klebsiella pneumoniae</i> (Carbapenem-resistant) <i>Pseudomonas aeruginosa</i> (Multidrug-resistant) <i>Staphylococcus aureus</i> (Methicillin-resistant) <i>Staphylococcus aureus</i> (Methicillin-resistant) <i>Staphylococcus epidermidis</i> (Methicillin-resistant)	2 minutes	Clean
EN 13727	Bactericidal	<i>Pseudomonas aeruginosa</i> <i>Staphylococcus aureus</i> <i>Enterococcus hirae</i> <i>Escherichia coli</i> <i>Enterococcus faecium</i> <i>Acinetobacter baumannii</i> <i>Klebsiella pneumoniae</i>	30 seconds	Clean
EN 13727	Bactericidal	<i>Staphylococcus aureus</i> <i>Escherichia coli</i> <i>Enterococcus faecalis</i> <i>Acinetobacter baumannii</i> <i>Klebsiella pneumoniae</i>	30 seconds	Dirty
EN 13727	Bactericidal	<i>Klebsiella pneumoniae</i>	1 minute	Dirty
EN 13727	Bactericidal	<i>Escherichia coli</i>	1 minute	Dirty
EN 13727	Bactericidal	Methicillin resistant <i>Staphylococcus aureus</i> Vancomycin resistant <i>Enterococcus faecalis</i>	1 minute	Clean
EN 13727	Bactericidal	<i>Acinetobacter baumannii</i>	1 minute	Dirty

Use biocides safely. Always read the product label and information before use.

Products within this range



Flowpack - 360 Wipes

Product code - S525110MPCE

Classification - Medical Device Class IIa

Packs per case - 6



Flowpack - 50 Wipes

Product code - S520110MPCE/S524110MPCE

Classification - Medical Device Class IIa

Packs per case - 24



Flowpack - 100 Wipes

Product code - S523110MPCE

Classification - Medical Device Class IIa

Packs per case - 12



Flowpack - 200 Wipes

Product code - S522110MPCE/S529110MPCE/ S625110MPCE

Classification - Medical Device Class IIa

Packs per case - 6



Canister - 200 Wipes

Product code - W615110MPCE

Classification - Medical Device Class IIa

Packs per case - 10



Canister - 240 Wipes

Product code - W520110MPCE

Classification - Medical Device Class IIa

Packs per case - 10

Products within this range



Bucket - 225 Wipes

Product code - W527110MPCE

Classification - Medical Device Class IIa

Packs per case - 4



Bucket - 350 Wipes

Product code - W522110MPCE

Classification - Medical Device Class IIa

Packs per case - 4



Refill - 225 Wipes

Product code - Q527110MPCE

Classification - Medical Device Class IIa

Packs per case - 4

Manufacturer:
Pal International Limited

Unit 3 Mountpark, Bardon II,
Victoria Road,
Ellistown,
Coalville,
LE67 1FA
United Kingdom

T: +44 (0) 1455 555700

E: sales@palinternational.com

W: www.palinternational.com

EU Representative:
Advena Ltd.,

Tower Business Centre
2nd Flr.
Tower Street
Swatar, BKR 4013
Malta

E: info@advena.mt

